ArtVentive Medical Group, Inc.
Endoluminal Occlusion System – EOS™ for Peripheral Use
Instructions for Use (IFU)

Caution: Federal (US) law restricts this device to sale by or on the order of a physician.

Device Name

ArtVentive (AVMG) Endoluminal Occlusion System – EOS™

Contents

- Delivery Catheter with preloaded implant consisting of a Nitinol scaffold and an ePTFE membrane
- 1 Patient Identification Card
- Guide Catheter and Dilator (Packaged Separately)

How Supplied

The ArtVentive Endoluminal Occlusion System - EOS™ is packaged and supplied sterile via EO sterilization and is non-pyrogenic. The Delivery Catheter with a preloaded Nitinol implant is packaged separately from the Guide Catheter. The System is intended for single use only.

Storage

Avoid exposure to water, sunlight, extreme temperatures and high humidity during storage. Store the EOS under controlled room temperatures. See the product label for the device shelf life. Do not use the device beyond the labeled shelf life.

Required Accessories Not Provided

- 1) 0.035" (or smaller) guide wire
- 2) 6Fr or larger Introducer sheath (for 3 - 8mm devices)
- 3) 8Fr or larger Introducer sheath (for 7.5 - 11mm devices)
- 4) Rotating hemostatic valves
- 5) Flushing syringe

Indication for Use

The ArtVentive Endoluminal Occlusion System - EOS™ is indicated for arterial and venous embolization in the peripheral vasculature.

Device Description

The ArtVentive Endoluminal Occlusion System - EOS™ device has been developed for percutaneous, peripheral occlusion of the peripheral arterial and venous vasculature. The system is comprised of the Delivery Catheter with a preloaded Implant and a Guide Catheter with Dilator (packaged separately).

The implant is made of a Nitinol scaffold with an impermeable ePTFE membrane and is designed with radial force sufficient to provide radial stiffness and strength against the vessel wall, and to minimize post-deployment migration. The EOS™ implant is available in three sizes serving vessel IDs from 3 - 11mm to accommodate the human vasculature (see Figure 1 for reference vessel diameters).

The EOS™ consists of the following components:

1) EOS™ Guide Catheter and Dilator

The Guide Catheter is a braided shaft with a stiff proximal section and a more flexible distal section to enable tracking through tortuous peripheral vasculature. A radiopaque marker on the distal end of the Guide Catheter is visible under fluoroscopy. The tip of the Guide Catheter is tapered to fit over the Dilator. The Dilator fits inside the Guide Catheter, exiting out through the distal end. The Dilator also has a tapered end for ease of advancement into the blood vessel. The guide wire and Dilator are removed from the Guide Catheter once it is in position for delivery of the implant.

a. Dimensions – Guide Catheter

- i. Total length: 92cm
- ii. Working length: 89cm

b. Dimensions – Dilator

- i. Total length: 108cm
- ii. Working length: 106cm

2) EOS™ Delivery Catheter

The EOS™ Delivery Catheter is compatible with the EOS™ Guide Catheter, and contains one (1) implant loaded on the distal end and a two-stage deployment handle on the proximal end. The catheter has a stiff proximal section for pushability and a flexible distal section for trackability. The deployment handle has a side port to accommodate syringe attachment to flush the catheter of air and to pre-expand the ePTFE membrane before deploying the implant. The handle has one (1) yellow and one (1) blue clip which prevent premature deployment of the implant. The handle has two (2) sliders for releasing the implant. The proximal slider releases the proximal end of the implant. The distal slider releases the distal end of the implant. The total length of the EOS™ Delivery Catheter is 117 cm with a working length of 105 cm.

3) EOS™ Implant

The EOS™ implant is a Nitinol scaffold with an ePTFE membrane to occlude blood flow. It is designed for sufficient radial strength against the vessel wall under normal blood pressure to minimize post-deployment migration. The EOS™ implant is available in three sizes:


Figure 1. Device Dimensions-

Contraindications for Use

- 1) Implantation in proximity with high motility (locomotor) muscles.

Warnings

- The safety and effectiveness of this device for cardiac uses (for example, patent ductus arteriosus or paravalvular leak closures) and neurological uses have not been established.
- This device should only be used by physicians who have been trained in transcatheter techniques. The physician should determine which patients are suitable candidates for procedures using this device.
- Do not use a power injection syringe to inject contrast solution through this device.
- Physicians must be prepared to deal with urgent situations, such as device embolization which require removal of the device. This includes the availability of an on-site surgeon.
- Do not use the device if the packaging sterile barrier is open or damaged.
- Use on or before the last day of the expiration month noted on the product packaging.
- The device is sterilized using ethylene oxide and is for single use only. Do not reuse or resterilize. Attempts to resterilize the device may result in device malfunction, inadequate sterilization, or patient harm.
- In the event that any resistance is felt during or following removal of the protective cover or prior to deployment, do not attempt to introduce the device. Premature deployment of the device could occur within the guide catheter prior to reaching the targeted treatment location. Select a new device.

Precautions

- Patients who are allergic to nickel may have an allergic reaction to this device.
- The physician should exercise clinical judgment in situations that involve the use of anticoagulants or antiplatelet drugs before, during and/or after use of the device.
- Use in specific populations
 - Pregnancy – care should be taken to minimize the radiation exposure to the fetus and the mother.

- Nursing mothers – there has been no quantitative assessment of the presence of leachables in breast milk.
- Only place device into vessel size as indicated for use. If placed in a vessel smaller than indicated for use, device may damage vessel wall. If placed in a vessel larger than indicated for use, device may migrate.

MRI Safety Information

MR Conditional

Non-clinical testing has demonstrated that the EOS Implant is MR Conditional for single lengths up to 20mm. A patient with this device can be scanned safely, immediately after placement under the following conditions:

- Static magnetic field of 1.5-Tesla and 3-Tesla, only
- Maximum spatial gradient magnetic field of 4,000-Gauss/cm (40 T/m) or less
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 4-W/kg for 15 minutes of scanning (i.e., per pulse sequence)
- First Level Controlled Mode of operation for the MR system

MRI-Related Heating

In non-clinical testing, the EOS implant produced the following temperature rises during MRI performed for 15-min of scanning (i.e., per pulse sequence) in 1.5-Tesla/64-MHz (Magnetom, Siemens Medical Solutions, Malvern, PA. Software Numaris/4, Version Syngo MR 2002B DHHS Active-shielded, horizontal field scanner) and 3-Tesla (3-Tesla/128-MHz, Excite, HDx, Software 14X.M5, General Electric Healthcare, Milwaukee, WI) MR systems:

	1.5-T, 8-mm	1.5-T, 5-mm	3-T, 8-mm	3-T, 5-mm
MR system reported, whole body averaged SAR	2.9-W/kg	2.9-W/kg	2.9-W/kg	2.9-W/kg
Calorimetry measured values, whole body averaged SAR	2.1-W/kg	2.1-W/kg	2.7-W/kg	2.7-W/kg
Highest temperature change	1.6°C	1.6°C	1.8°C	1.7°C
Temperature scaled to whole body averaged SAR of 4-W/kg	2.2°C	2.2°C	2.5°C	2.3°C

Artifact Information

The maximum artifact size as seen on the gradient echo pulse sequence at 3-Tesla extends approximately 5-mm relative to the size and shape of this implant when imaged with a gradient echo pulse sequence and a 3 Tesla MRI system. The artifact does obscure the device lumen.

Procedural Steps

Preparation

- 1) Obtain access to the vasculature using standard transcatheter techniques.
- 2) Perform an angiogram to measure the vessel diameter at the desired occlusion site.
- 3) Inspect the device for kinks, breaks or for premature deployment of the scaffold. Inspect sterile packaging to verify that no damage has occurred as a result of shipping. Do not use if the device or the sterilization barrier has been damaged or broken.
- 4) Remove the Guide Catheter from the package and remove the dilator from the Guide Catheter.
- 5) Flush all components of the system to remove air.

Advancing the EOS™ Guide Catheter

- 1) Advance Guide Catheter with dilator over a .035" or smaller guide wire through the introducer sheath.
- 2) Using the radiopaque marker as a guide, position the Guide Catheter at the targeted implantation site.
- 3) Remove the guide wire and the dilator while keeping the Guide Catheter in place.

Delivering the EOS™ Implant

- 1) Remove the Delivery Catheter with implant from the packaging hoop.
CAUTION: Inspect the implant for any sign of damage or partial deployment and replace with a new device if necessary. Do not attempt to replace the protective cover if the implant has partially deployed.
- 2) Attach a 10cc luer-lock syringe of saline to the flush port on the deployment handle.

- 3) With the protective cover on the EOS™ device, flush the delivery system to remove air.
- 4) With the protective cover in place, position the EOS™ device at the hub of the Guide Catheter.
- 5) Advance the EOS™ device through the protective cover into the Guide Catheter. Once the EOS™ device is completely in the Guide Catheter, slide the protective cover proximally over the delivery catheter toward the deployment handle.
- 6) Advance the EOS™ delivery system to the desired treatment location.
- 7) Pull the Guide Catheter back to expose the EOS™ device. Confirm with fluoroscopy that the EOS™ device is completely beyond the tip of the Guide Catheter and positioned for deployment.
- 8) Remove the yellow proximal clip from the handle.
- 9) Attach 10cc luer-lock syringe of 80% saline and 20% contrast to the flush port. Under fluoroscopy, flush the device with a powerful flush to fully expand the ePTFE membrane prior to deployment. It is critical at this point to maintain a positive pressure and fully expand the membrane with continued flush.
- 10) To deploy the proximal end of the EOS™ device, continue flushing while pulling the proximal slider. It is important to continue to flush while pulling the proximal slider. NOTE: Once the proximal end of the EOS™ device has been deployed, the device can no longer be recaptured and the full delivery sequence must be completed.
- 11) Evaluate the device position with fluoroscopy. Note: The device may be repositioned at the target site at this point, if desired. NOTE: The EOS™ device may be repositioned distally at the target site at this point, if desired. To reposition, gently advance the delivery catheter to move the EOS™ device distally. Do not attempt to retract the EOS™ device, as this may damage the device and prevent it from deploying correctly.
- 12) Remove the distal blue clip from the deployment handle.
- 13) To complete deployment, continue flushing while pulling the distal slider to deploy the distal end of the EOS™ device. It is important to continue to flush while pulling the distal slider. The implant is now fully deployed and detached from the Delivery Catheter.
- 14) Evaluate device placement and blood flow under fluoroscopy for desired results.

Symbols

Symbol	Definition
	Do not Reuse - Single use only
	Use by expiration date on the package
	Method of sterilization using Ethylene Oxide
	Consult instructions for use
	Manufacturer
	Do not use the device if the packaging has been damaged or if sterility has been compromised
	MR Conditional
	Do not resterilize
	Keep away from sunlight
	Keep dry
	Non-pyrogenic

Disposal

Dispose of the Guide Catheter, Dilator, and Delivery Catheter per institutional standard for biohazard materials.

Post Procedural Instructions

1. Patient ID Card – A patient identification card is enclosed in the package. Complete this card and give to the patient.
2. MedicAlert Service – Recommend the patient to become a MedicAlert member by calling 1.888.633.4298 or enrolling online at www.medicalert.org. The patient will be asked to provide his/her condition, implanted device name, and restrictions concerning use of MRI.

Return of Devices

If any portion of the ArtVentive Endoluminal Occlusion System™ fails prior to or during a procedure, discontinue use. The device should be placed in a biohazard container with a Returned Material Authorization (RMA) number assigned by ArtVentive. Please call your local distributor or AVMG Customer Service at +1 650 465-5259 for an RMA number.

ArtVentive Medical Group, Inc.
1797 Playa Vista
San Marcos, California 92078
USA
Tel: +1 (650) 465-5259
Fax: +1 (650) 240-1706
customerservice@artventivemedical.com

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Warranty

ArtVentive Medical, Inc. warrants that reasonable care has been used in the design and manufacture of this device. This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether expressed or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness for a particular use. Handling, storage, cleaning, and sterilization of this device as well as other factors relating to the patient, diagnosis, treatment, surgical procedures, and other matters beyond ArtVentive's control directly affect the device and the results obtained from its use. ArtVentive's obligation under this warranty is limited to the repair or replacement of this device and ArtVentive shall not be liable for any incidental or consequential loss, damage, or expense directly or indirectly arising from the use of this device. ArtVentive neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this device. ArtVentive assumes no liability with respect to devices reused, reprocessed, or resterilized and makes no warranties, expressed or implied, including but not limited to merchantability or fitness for intended use, with respect to such device.

ArtVentive Medical Group, Inc.
Endoluminal Occlusion System – EOS™ para uso periférico
Instrucciones de uso (IFU)

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ArtVentre Medical Group, Inc.
MDSS GmbH