# MAZIC<sup>®</sup> Duro

Nano Hybrid Ceramic for CAD/CAM

 Instruction for Use

## Description

MAZIC<sup>®</sup> Duro, a dental restorative material, is a radiopaque, block type and disk, hybrid ceramic that is intended for the fabrication of inlays, onlays, veneers and crowns. It can be matched to the patient's tooth shade and is intended for use in dental CAD/CAM milling procedures.

## Composition

The polymeric matrix consists of a high degree of cross linking monomer. The inorganic fillers are barium aluminosilicate, silicon dioxide and zirconia. The total content of inorganic fillers is approximately 80wt% and particle size of inorganic fillers is between 0.01  $\mu\text{m}$  and 1  $\mu\text{m}$ . Additives are initiators, stabilizers, pigments, etc.

## Indications

Fabrication of inlays, onlays, veneers and crowns

## Recommended procedure

1. Remove the inappropriate part of teeth/tooth for dental treatment, carefully.
2. If pulp exposure concerns, use a minimum amount of calcium hydroxide.
3. Place a temporary restoration on tooth before permanent restoration.
4. Select the proper shade of the mill block to be based on patient's teeth.
5. Fix the block / Disc at milling device(CAD/CAM) with an appropriate holder(jig).
6. Mill the block / Disc to proper size by using CAD/CAM
7. Remove temporary restoration and clean tooth completely.
8. As required, apply a suitable bonding system such as etchant, primer and adhesive to the prepared tooth and restoration according to manufacturer's recommendation.
9. Apply a resin cement according to manufacturer's recommendation and place(seat) a milled block on tooth.
10. Remove excess resin cement and finish with polishing.


## Storage

1. Avoid exposure to the sunlight.
2. Store at room temperature.

## Precautions

1. It should be used by a dentist and dental professional. Observe the procedure and the precautions.
2. Keep it away from a reach of non-dental user including children, senior citizens and others.
3. If do not keep the following recommendation, it may not get an ideal results.
  - Framework is not to be thin than required minimum thickness for restoration.
  - Mill the block by compatible CAD/CAM system.
  - Select the proper bur to this hybrid ceramic block.
  - Do not apply strong torsion and shear to the handle intentionally.
4. Do not breathe in dust particles. Use of air exhauster, mask, protective gloves and glasses is recommended.
5. Do not use after expiry date.
6. If skin sensitization or allergy(methacrylate, etc) occurs to this product, stop using.
7. Do not store or exposure at elevated temperature or intense light.
8. Apply for use purpose only.
9. Be careful not to fall down the product with mandrel.
10. If contaminant occurs to the product, do not use.
11. For a single use only.

## Symbols

## Warranty

VERICOM CO., LTD. will replace product that is proven to be defective. VERICOM CO., LTD. does not accept liability for any damage or loss, direct or consequential, stemming from the use of or inability to use product described. Before using, it is the responsibility of the user to determine the suitability of the product for its intended use. The user assumes all risk and liability in connection therewith.

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