


eazyplex[®] SARS-CoV-2


When it really counts — eazyplex[®] SARS-CoV-2!

The best combination of sensitivity, specificity and speed!




- Lyophilized tests
- 1 min sample pretreatment
- No RNA extraction
- Variable device selection
- Real-Time-Detection
- **CE IVD**

Take advantage of the optimal combination of the most innovative nucleic acid amplification technology and the user-friendly format of the eazyplex[®] system. With the eazyplex[®] SARS-CoV-2 you will receive the result without RNA extraction directly from UTM/VTM, eSwab medium or even dry swabs in 25 minutes at the latest. Our devices guide you through the test with their intuitive software. The sample control informs you about the swab quality and the automatically generated table provides clear results. Maximum diagnostic reliability is ensured by the dual testing of the WHO-recommended target "N gene" and the ORF8 gene.



Find out more and contact us at support@eazyplex.com



The **eazyplex[®]** test system is a freeze dried ready to use amplification reaction and can be stored direct at the workplace. During the amplification the result is shown in Real-Time.



AmplexDiagnostics GmbH
Werkstraße 2
83555 Gars Bahnhof
Germany
www.eazyplex.com

Kontakt/Contact

info@eazyplex.com
Tel.: +49 (0) 8073/91693-50
Fax: +49 (0) 8073/91693-55

