AYUDAS

813

Rubber Tip for canes.

Sizes: 10-12-14-16-18-20-22mm



814

Rubber Tips with suction cups. Size: 16 -18 mm.



815

Black Rubber Tip for Crutches. Metal ring in the inside.

Size: 19 - 21 mm.



816

Rubber Tip for canes. Sizes: 10-12-14-16-18-20-22mm.



817

Rubber Tips with suction cups. Size: 16 -18 mm.



818

Black Rubber Tip for Crutches. Metal ring in the inside. Size: 19 - 21 mm.



822

Rubber Tips for walker. Metal Ring inside Sizes: 22-25 y 28 mm




823

Rubber tips for shower bench. Metal ring inside. Sizes: 25 - 28 mm.



824

Rubber Tips for peddlers





Rubber Tip for aluminium cane. Metal ring inside. Size: 14mm



826

Brown rubber Tip for aluminium cane. Metal ring inside. Size: 14mm



827

Illuminated Tip,very useful to walk during the night or when it is quite dark. Metal ring inside. Size: 19 mm



828

Trekking rubber tip. Metal ring inside. Size: 12mm



834

Hexagonal rubber tip. Metal ring inside. Size: 21mm



8151191

Colourful rubbers for crutches. Colours: pink, purple, red and green. Size: 19mm.



831

Rubber tip with bellow. Base: 4 cm Size: 16 and 18 mm



Rubber tips with wide base for all types of canes (aluminum and wood). They provide greater stability and security in support. In addition, they allow the cane to remain standing.














820122 Black Beach tip 22 mm

82019 Grey Beach tip 19 mm

Grey beach tip 22 mm



Its great flexibility allows an important part of its base to always support. Its special shape achieves a high degree of cushioning, reducing the injuries of elbows and shoulders of its users.



819

Rubber with bellows. It has a great capacity of fl exibility that allows a constant contact of the base to the ground, achieving a high adhesion and great capacity of absorption to the impacts. Base diameter 5.5 cm.

Sizes: 17,19 and 21 mm



Rubber tip. Its design allows the cane to remain vertical in most situations. Its large base gives more security and stability preventing the cane from sinking in sand, mud or wet areas. Made of rubber and with an inner metal ring. Metal base inside.



8291

Tripod Rubber. Base: 7,5 cm

Sizes: 16 and 19 mm





829

Tripod Rubber. Base: 10 cm

Size: 16 and 19 mm



840 / 842 / 844 / 846

Its base allows the user to be safe in every type of ground or surface. The cane remains standing.

AYUDAS MOVILIDAD M. GARCIA

840

Quad tip base 8,6 x 8,6 cm. For aluminum canes. Sizes: 16mm and 19 mm.



846

Suspension rubber tip with three support points. Rotatory. It includes two plastic pieces to fit in 16mm and 19mm.



842

Quad tip base 15x15 cm. For aluminum canes. Size: 16mm and 19mm.



844

Cane tip Wide base for every type of canes (aluminum and wood). Sizes: 16mm and 19mm.





Cane plastic holder for walkers and wheelchairs.



837

Rotative plastic cane holder.





810

Black, brown, garnet or beige cord handle. Very easy use.




811

Black cover for crutches holder. Made of neoprene. Very soft.



830

Supporting cane clip. Made of nylon. This way, the cane will not fall. It will be supported on the table, chair, etc



Display for 10 folding canes. Made of acrylic.

795

Carton display for 12 units.



27 cm. 28 cm.



ACCESSORIES

Three levels wooden display for 15 pieces with wheels. Size:46x21x71cm

Round wooden display for 18 pieces. Size: 33 x 50 cm
Round wooden display for 12 pieces. Size: 27 x 50 cm





Steel rotatory display. 3 levels. 42 units



797 800 801

Steel Wire display for walking sticks and smaller pieces (folding canes and accesories). Available in black and white.



Dimensions: 40 x 183 cm.





Rotary display for folding canes. Diameter: 40 cm. Height: 41 cm.







Dimensions:

Width: 40 cm

Height: 183 cm



809

Dimensions: Width 60 cm Height: 183 cm

808 809

Steel Wire Display for smaller pieces (blister, accessories, etc.) Available in black and white