

# MAZIC D Material Compatibility Table\_V6



	Material										Material								
	Machine Name	TEMP	SG	Mode 2.0	Ortho Model	Cast	Splint	Denture	Etc.		Machine Name	TEMP	SG	Model 2.0	Ortho Model	Cast	Splint	Denture	Etc.
DLP	MAZIC D P500	●	●	●	●	●			Machine Setting		Sindoh A1+Plus	Laser Power :28 Default	Laser Power :32 Default	Laser Power :25 Default	Laser Power :25 Default				
	Hunter/p100	●	●	●	●	●			SLA	Formlab 2	Denture Teeth A2 50um								
	Carima 1M2	●	●	●	●	●				Zenith U	Dreve_Castable_SoftVAT	Clear_SG_SoftVAT	Dreve_Castable_SoftVAT		Castable Red_SoftVAT				
	Carima 1M3	●	●	●	●	●			Machine Setting	Karv1. 2		●	●	●	●	●			
	Vetz 30 D2	●	●	●	●	●				Ackuretta		●	●		●	●			
	Asiga-MAX UV	●	●	●	●	●				Freeshape 120		●	●		●	●			
	Zenith D	●	●			●				Phrozen		●	●	●					
	B103D W11	Nextdent SG & Ortho	Nextdent SG & Ortho	Nextdent Model B : 80s O : 8s	Nextdent SG & Ortho	Nextdent SG & Ortho				Soinc Mini		●	●	●					
	envisionTEC E1	●	●							Sonic 2K		●	●	●	●	●			
	Nextdent 5100	Nextdent MFH & Ortho Rigid	Nextdent SG	Nextdent Model 2.0 Gray	Nextdent Ortho Rigid					XL 4K		●	●	●	●	●			
	Kulzer cara 4.0	dima Impression Blue 100um	dima guide 100um	dima Print Model HD Standard 50um	dima guide 100um					Shuffle XL		●	●	●	●	●			
										Shuffle K. Edition			15-16s /100um	●					
										Shuffle 2K		●	15-16s /100um						
										Anycubic Photon Zero	About 20-25s /50um	About 10-12s /50um			About 10-12s /50um				
										Photon	●	●	●		●				
									Photon S	●	●	●		●					
									Photon Mono X	50um B : 40s O : 7s	50um B : 35s O : 6s								
									Uniz Slash 2	●	●								
									Shash+		●								
									Dazz L120 Zeron3D Pixel one	●	●	●							
									WANHAO Duplicator D7		●	●							



Available



In Progress



Double click & Download



Closed Type